



stepahead[®]

Comfort

Stepahead Comfort is designed to cushion and protect the foot. Ideal for the older foot type where the fibro-fatty pad is thin or where plantar lesions such as corns and calluses need pressure deflection.

Image for illustration purposes only.

When to specify?

The compressible properties of EVA make it an ideal choice when a more progressive type of orthotic control is required. Rigid or immobile feet often prefer this approach and patient compliance is often improved.

When treating?

- Foot conditions associated with excessive pronation or supination
- Excessive 'shock' or loading during gait

In which shoes?

- Most casual styles*
- Trainers*
- Walking boots*

*EVA devices tend to take more room in the shoes than polypropylene devices and consequently this must be taken into consideration when prescribing. Please ask for advice from our technical team.



techstep

TECHNOLOGY THAT'S A STEP AHEAD



Comfort

RM763530
10/07/0.5K