



stepahead®

Dynamic

Stepahead Dynamic is aimed at the more active foot where a degree of impact attenuation is needed along with functional control at the propulsive phase of gait i.e. under the metatarsal heads and toes.

Image for illustration purposes only.

When to specify?

Where the foot/orthotic interface needs to be applied to the whole of the foot – heel to toes. The 4 different densities of EVA available allow the support offered to be specific to the individual with multi-density devices also being available.

When treating?

- Foot conditions that require support and control and all stages of the gait cycle
- Active foot types with propulsive phase control needs, eg. runners

In which shoes?

- Most casual styles*
- Trainers and walking boots*

*EVA devices tend to take more room in the shoes than polypropylene devices and consequently this must be taken into consideration when prescribing. Please ask for advice from our technical team.



SALTS
techstep

TECHNOLOGY THAT'S A STEP AHEAD



Dynamic

RM763531
10/07/0.5K