



stepahead®

## Sprint

Stepahead Sprint is specifically designed for propulsive phase control as needed by cyclists and sprinters. It allows functional control to be prescribed and integrated into the orthotic to effect this phase of the gait cycle and at the same time fit the footwear i.e. running spikes and cycle shoes.

Image for illustration purposes only.

### When to specify?

Where propulsive phase control is required within the confines of sports specific footwear that has little room for conventional orthotic devices.

### When treating?

- Pathologies associated with metatarsalphalangeal dysfunction during the propulsive phase of gait

### In which shoes?

- Cycling cleats
- Football and rugby boots
- Running spikes
- Trainers



SALTS  
techstep

TECHNOLOGY THAT'S A STEP AHEAD



Sprint

RM763534  
10/07/0.5K